




Accessories & Spares

for Rocket Synergy Evaporators






Disclaimer

The original, authoritative version of this document is the English version produced by Genevac Ltd. Subsequent changes to any document made by any third party have not been reviewed nor authenticated by Genevac Ltd. Such changes may include, but are not limited to, translation into languages other than English and additions to or deletions from the original content. Genevac Ltd disclaims responsibility for any and all damages, injuries, warranty claims and/or any other claims associated with such changes, inasmuch as such changes result in content that is different from the authoritative Genevac published version of the document.

Table of Contents

Introduction	4
Limitations of Use	4
Surface Finish Variation	4
Sample Holder Maintenance	4
Glassware	5
Labelling of Glassware	5
Seal Types	6
Acceptable Solvents	7
Rotors	8
Rotor Spares	9
Rotor Accessories	11
Sample Holders for Flask Rotors	12
Standard Flask	12
Standard Flask Spares	13
GC Vial System ²	14
GC Vial System ² Spares	15
SampleGenie ⁴	18
ASE [®] Tube Puck	24
Flip-Flop	26
Dionex™ 250 ml ASE [®] Bottle Adaptor	29
Rocket Synergy - Accessories	30
Infinity Trolley	31
Rocket Synergy - General Spares	32
Pre-Synergy Rockets - Rotor Compatibility	34



Introduction

Genevac produce a range of sample holders, rotors and accessories that can be used to enhance performance and increase the flexibility of Rocket Synergy and Rocket 4D Synergy evaporators.

Sample holders are precision manufactured using high-quality materials such as borosilicate glass and high-grade aluminium to ensure the best possible evaporation.

Limitations of Use

This brochure contains information on a range of accessories and sample holders intended for use with Rocket Synergy and Rocket 4D Synergy evaporators. Some accessories and sample holders are compatible with earlier versions of the Rocket evaporator, although limitations of use may apply. Where limitations of use do apply, they will be indicated as shown in the following example:

For use with Rocket 4D Synergy only.

In this example, the accessory must not be used with any system other than that specified.

Some sample holder may not be compatible with older model Rocket evaporators. Refer to ***Pre Synergy Rockets - Rotor Compatibility*** in this brochure.

Surface Finish Variation

For some components, the colouring of the surface finish may vary slightly from batch to batch - for example, one batch of vial adaptors may appear slightly lighter or darker than another.

Variation in the colour of the surface finish will not affect the performance or integrity of the component, or the assembly.

Sample Holder Maintenance

The Rocket evaporator controls pressure in the inner and outer vacuum chambers independently. Sample holders form a vital part of the seal between these chambers, therefore it is important that a perfect seal be formed around each sample holder. To achieve this, it is critical that sample holders are maintained in good condition and are used correctly.

Refer to the ***Genevac Rocket Sample Holders User Manual*** for advice on day-to-day maintenance and use of Rocket sample holders.

Glassware

Some Genevac evaporators subject samples to accelerations in excess of 500 times gravity. This ensures excellent evaporation, but places a high degree of stress on the glassware. Selection of suitable glassware is essential to avoid glassware breakage.

Always use high-quality glassware from recognised manufacturers of international repute. Experience has shown that 'own-brand' catalogue glassware can have highly variable dimensions and manufacturing quality; this kind of glassware is not recommended for use in Genevac evaporators due to excessive risk of breakage.

When selecting glassware, ensure the following:

- Reputable scientific glassware manufacturer
- High quality, machine-made borosilicate glass
- Minimum 1.2 mm wall thickness

Re-use of glassware in centrifugal evaporators gradually reduces the strength of the glassware. Replace all glassware regularly to avoid breakage.

Damaged Glassware

Even a small scratch can reduce the strength of glassware by a significant degree. **Do not use glassware that is scratched, chipped, cracked or otherwise damaged. Do not use metal spatulas, wire brushes or similar hard instruments to clean glassware or remove dried samples.**

Labelling of Glassware

There are two common methods of labelling glassware in laboratory environments. Both types can cause problems in the evaporator.

Paper/Thin Film Labels

The label effectively increases the diameter of the tube or flask (especially if it wraps around more than once). When the tube or flask is not fully snug in the holder, heat transfer is reduced and any raised points (e.g. tiny fragments of glass from a previous breakage, small lumps of paper towel left behind after cleaning, or the label itself) place additional stress on the glassware.

Ceramic Etched

Ceramic-etched labelling weakens the glass. Genevac recommend testing new tubes for breakage before purchasing glassware in bulk.

Seal Types

Seals for a selection of Rocket Synergy sample holders are available in two variants, to suit different solvent applications.

Viton®

For less aggressive applications, low-leachable white Viton seals are available.

Perfluoroelastomer

For applications in which Viton seals are unsuitable, highly chemical resistant perfluoroelastomer seals are available.


Compatibility

To ensure compatibility with your application, please refer to the **Acceptable Solvents** table on the following page before purchasing. If your solvent is not listed, please contact your local Genevac representative.


Identifying Perfluoroelastomer Seals

One-piece perfluoroelastomer seals (GC Vial System², SampleGenie⁴, Flip-Flop) are identified by a small notch in the periphery of the seal.

Two-piece vial seals (SampleGenie⁴ only) consist of a holder and a split insert. The perfluoroelastomer version is indicated by a green mark on the periphery of the holder.



One-piece perfluoroelastomer vial seal showing indicator notch.



Two-piece perfluoroelastomer vial seal showing green indicator mark.


Acceptable Solvents


Solvent	Standard Flask	GC Vial System ²	SampleGenie ⁴	ASE Tube Puck	Flip-Flop	ASE Bottle Adaptor
Acetic Acid	■	□	□	■	■	■
Acetic Anhydride	■	□	□	■	□	■
Acetonitrile	■	□	□	■	□	■
Ammonium Hydroxide	■	□	□	■	□	■
1-Butanol	■	■	■	■	■	■
2-Butanol	■	■	■	■	■	■
2-Butanone	■	□	□	■	□	■
Butyl Acetate	■	■	■	■	■	■
<i>tert</i> -Butyl Alcohol	■	■	■	■	■	■
Chlorobenzene	■	■	■	■	■	■
Chloroform	■	■	■	■	■	■
Cyclohexane	■	■	■	■	■	■
1,2-Dichloroethane	■	□	□	■	□	■
Dichloromethane	■	□	□	■	□	■
Diethyl ether	■	□	□	■	□	■
Diisopropyl ether	■	□	□	■	□	■
N,N-Dimethylformamide	■	□	□	■	□	■
1,4-Dioxane	■	□	□	■	□	■
Ethanol	■	□	□	■	□	■
Ethyl Acetate	■	□	□	■	□	■
Formic Acid	■	□	□	■	□	■
Heptane	■	■	■	■	■	■
Hexane	■	■	■	■	■	■
Methanol	■	□	□	■	□	■
Methyl- <i>tert</i> -butyl ether	■	□	□	■	□	■
Pentane	■	□	□	■	□	■
1-Pentanol	■	□	□	■	□	■
Petroleum Ether	■	■	■	■	■	■
1-Propanol	■	■	■	■	■	■
2-Propanol	■	■	■	■	■	■
Propanone	■	□	□	■	□	■
Pyridine	■	□	□	■	□	■
Sodium hydroxide (aqueous)	■	□	□	■	□	■
Tetrahydrofuran	■	□	□	■	□	■
Toluene	■	□	□	■	□	■
Triethylamine	■	□	□	■	□	■
Trifluoroacetic acid	■	□	□	■	□	■
2,2,2-Trifluoroethanol	■	□	□	■	□	■
Water	■	■	■	■	■	■
Common Mixtures						
Water & Acetonitrile	■	□	□	■	□	■
Water & Methanol	■	□	□	■	□	■
Water & Ammonia	■	□	□	■	□	■
Methanol & DCM	■	□	□	■	□	■
TFA & DCM	■	□	□	■	□	■


- Solvent may be evaporated using sample holder.
- Only acceptable if sample holder has perfluoroelastomer seals.

NOTE: Consult your local Genevac representative for advice before using sample holders with solvents that are not listed above.


Rotors


	Description	Flask Rotor	Part Number 70-1734/S
	Details	Six-position Flask Rotor assembly for Rocket Synergy and Rocket 4D Synergy evaporators. Includes inner lid and six blanking plugs.	


	Description	Bowl Rotor - Batch	Part Number 70-1820/S
	Details	Five-litre bowl rotor for Rocket Synergy or Rocket 4D Synergy. Includes batch inner lid (no feed), batch baffles, lifting handles and USB data-key with Rocket Synergy Batch Methods.	


	Description	Bowl Rotor - Feed	Part Number 70-2041/S
	Details	Five-litre bowl rotor for Rocket 4D Synergy (only). Includes feed inner lid, feed baffle, lifting handles and USB data-key with Rocket Synergy Feed Methods.	


Rotor Spares


	Description	Inner Lid - Non Feed	Part Number
	Details	Inner lid for Flask Rotor or Batch (non-feed) Bowl Rotor.	
			70-1548/S


	Description	Inner Lid - Feed	Part Number
	Details	Inner lid for feed Bowl Rotor. For use with Rocket 4D Synergy feed systems only.	
			70-1547/S


	Description	Batch Baffles	Part Number
	Details	Set of two batch baffles for Bowl Rotor.	
			70-1549/S


	Description	Feed Baffle	Part Number
	Details	Feed baffle for Bowl Rotor.	
			04-7652/S

	Description	Inner Lid Seal	Part Number 04-7448/S
	Details	Inner lid seal for Flask & Bowl Rotors.	


	Description	Flask Seals	Part Number 70-1180
	Details	Set of six Flask Seals. Compatible with all Rocket Synergy Flask Rotors.	


	Description	Rotor Hub O-Ring	Part Number 04-5676
	Details	O-ring (Viton) for Rotor Hub.	


	Description	Pouring Spout Cap Seal	Part Number 70-1556/S
	Details	Seal kit (consisting of plug & O-ring) for Bowl Rotor pouring spout cap.	


	Description	Pouring Spout Cap	Part Number 70-1557/S
	Details	Cap for Bowl Rotor pouring spout.	


Rotor Accessories

	Description	Lid Stand	Part Number
	Details	Stand for Inner lid. Suitable for use with all variants.	
			70-1186

	Description	Lifting Handles	Part Number
	Details	Set of two lifting handles for Rocket Synergy Bowl Rotor.	
			70-1552/S

	Description	Blanking Plug	Part Number
	Details	Set of two blanking plugs for Flask Rotor.	
			70-2004/S


	Description	Pouring Stand	Part Number
	Details	Pouring Stand. Suitable for use with all Bowl Rotors with pouring spout.	
			70-1798

	Description	Synergy Methods	Part Number
	Details	USB data-key containing evaporation method files for Rocket Synergy Flask Rotor and Bowl Rotor (feed and/or batch) methods.	
			04-8844


Sample Holders for Flask Rotors


Standard Flask




	Description	Standard Flask Assembly	Part Number
	Details	Standard Flask with collar pad and collar. Set of six.	
			70-1321


Standard Flask Spares


	Description	Standard Flask	Part Number
			70-1170
	Details	Standard Flask - flask only. Set of six.	


	Description	Collar Pad	Part Number
			70-1317
	Details	Collar Pad for Standard Flask. Set of six. <i>For evaporators with serial numbers lower than RK0136, flask adaptors must be removed from the rotor.</i>	

	Description	Collar	Part Number
			70-1319
	Details	Collar for standard flask. Set of six. <i>For evaporators with serial numbers lower than RK0136, flask adaptors must be removed from the rotor.</i>	


GC Vial System²





	Description	GC Vial System ² (Viton Seals)	Part Number 70-1423
	Details	GC Vial System ² flask assembly with white Viton seals. Set of six.	


	Description	GC Vial System ² (Perfluoroelastomer Seals)	Part Number 70-1424
	Details	GC Vial System ² flask assembly with perfluoroelastomer seals. Set of six.	


GC Vial System² Spares


	Description	GC Vial System ² Flask	Part Number
	Details	Flask only, for GC Vial System ² . Set of six.	
			70-1425


	Description	Vial Seal (White Viton)	Part Number
	Details	White Viton vial seal for GC Vial System ² . Set of six.	
			70-1177


	Description	Vial Seal (Perfluoroelastomer)	Part Number
	Details	Perfluoroelastomer vial seal for GC Vial System ² . Set of six. NOTE: Notch in periphery indicates perfluoroelastomer.	
			70-1178

	Description	GC Vial Adaptor ²	Part Number
	Details	Two-piece vial adaptor with O-ring for GC Vial System ² . Set of six. NOTE: Backwards compatible with all GC Vial System ² flask assemblies. Replaces one-piece vial adaptor (70-1428).	
			70-1800

	Description	Vial Adaptor O-ring	Part Number
	Details	O-ring for GC Vial System ² Vial Adaptor (70-1800). Set of six.	
			70-1801

	Description	GC Vial System ² Collar Pad	Part Number
	Details	Collar Pad for GC Vial System ² . Set of six.	
			70-1317


	Description	GC Vial System ² Collar	Part Number
	Details	Collar for GC Vial System ² . Set of six.	
			70-1320


	Description	GC Vial System ² Cap	Part Number
	Details	Cap for GC Vial System ² . Set of six.	
			70-1176

This page deliberately left blank.




SampleGenie⁴




	Description	SampleGenie ⁴	Part Number
	Details	70-1740 SampleGenie ⁴ Flask assembly. Set of six. For use in Rocket Synergy evaporators and Rocket 4D Synergy evaporators with optional Flask Rotor. Evaporators with serial numbers lower than RK0136, flask adaptors must be removed from the rotor.	
NOTE: Excludes vial adaptor, vial adaptor pad and vial seal. When ordering, please select vial adaptor and vial seals from the tables on the following pages.			


SampleGenie⁴ - Vial Adaptors


	Description	Vial Adaptor	Part Number
	Details	Vial Adaptor for SampleGenie ⁴ . Set of six. NOTE: Includes Vial Adaptor Pad.	


	Max Vial Diameter (mm)	Vial Length (mm)	Seal Required?	Part Number
DG	15.3	38.5 to 41.5	One-part	70-1766
DD	15.3	43.5 to 56.5	One-part	70-1763
DK	16.5	46.5 to 51.5	One-part	70-2028
DJ	17.3	45.5 to 51.5	Two-part	70-1769
DC	17.3	58.5 to 61.5	One-part	70-1762
DH	20.8	41.5 to 44.5	Two-part	70-1767
DF	23.3	44.5 to 47.5	Two-part	70-1765
DE	23.3	66.5 to 71.5	Two-part	70-1764
DI	23.8	56.5 to 59.5	Two-part	70-1767
DB	28.0	58.5 to 61.5	Two-part	70-1761
DA	28.0	70 (max)	Two-part	70-1760


SampleGenie⁴ - Vial Seals


	Description	Two-Part Vial Seal - Insert	Part Number 70-1757
	Details	Seal Insert for SampleGenie ⁴ two-part vial seal. Set of six. For use with Viton vial seal (70-1758) or perfluoroelastomer vial seal (70-1759).	


	Description	Two-Part Vial Seal	Part Number 70-1758
	Details	Viton Vial Seal for SampleGenie ⁴ two-part vial seal. Set of six. Requires Seal Insert (70-1757).	


	Description	Two-Part Vial Seal	Part Number 70-1759
	Details	Perfluoroelastomer Vial Seal for SampleGenie ⁴ two-part vial seal. Set of six. Requires Seal Insert (70-1757). NOTE: Green mark in periphery indicates perfluoroelastomer.	


	Description	One-Part Vial Seal (White Viton)	Part Number 70-1196
	Details	One-part vial seal for SampleGenie ⁴ in white Viton.	

	Description	One-Part Vial Seal (Perfluoroelastomer)	Part Number 70-1197
	Details	One-part vial seal for SampleGenie ⁴ in perfluoroelastomer. NOTE: Notch in periphery indicates perfluoroelastomer.	


	Description	Vial Adaptor Pad 16 mm	Part Number 70-1745
	Details	16 mm Vial Adaptor Pad for SampleGenie ⁴ . Set of six.	


	Description	Vial Adaptor Pad 18 mm	Part Number 70-1746
	Details	18 mm Vial Adaptor Pad for SampleGenie ⁴ . Set of six.	


	Description	Vial Adaptor Pad 24 mm	Part Number 70-1747
	Details	24 mm Vial Adaptor Pad for SampleGenie ⁴ . Set of six.	


	Description	Vial Adaptor Pad 28 mm	Part Number 70-1748
	Details	28 mm Vial Adaptor Pad for SampleGenie ⁴ . Set of six.	

SampleGenie⁴ Spares

	Description	SampleGenie ⁴ Flask	Part Number
			70-1741
	Details	Flask assembly (only) for SampleGenie ⁴ . Set of six.	

	Description	Collar Pad	Part Number
			70-1743
	Details	Collar Pad (only) for SampleGenie ⁴ . Set of six.	

	Description	Collar	Part Number
			70-1742
	Details	Collar (only) for SampleGenie ⁴ . Set of six.	


	Description	Holder	Part Number
			70-1744
	Details	Holder for SampleGenie ⁴ Vial Adaptor. Set of six.	

This page deliberately left blank.




ASE® Tube Puck





	Description	ASE Tube Puck	Part Number
	Details	ASE Tube Puck assembly. Allows up to 18 Dionex ASE tubes to be processed simultaneously. Set of six.	

ASE Tube Puck Spares

	Description	Puck Seal	Part Number
	Details	Seal (Viton O-ring) for ASE Tube Puck. Set of 18.	
			70-1287


40 ml ASE Tube Adaptors


	Description	40 ml ASE Tube Adaptor (White Viton)	Part Number
	Details	40 ml ASE Tube Adaptor with white Viton seals. For use with ASE Tube Puck. Set of 18.	
			70-2008

	Description	40 ml ASE Tube Adaptor (Perfluoroelastomer)	Part Number
	Details	40 ml ASE Tube Adaptor with perfluoroelastomer seals. For use with ASE Tube Puck. Set of 18.	
			70-2009


Flip-Flop





	Description	Flip-Flop (Viton Seals)	Part Number
			70-2035
	Details	Flip-Flop assembly with white Viton seals. For use with ASE Tube Puck. Set of 18.	


	Description	ASE Tube Puck	Part Number
			70-2036
	Details	Flip-Flop assembly with perfluoroelastomer seals. For use with ASE Tube Puck. Set of 18.	


Flip-Flop Spares


	Description	Flip-Flop Tube	Part Number
			70-1297
	Details	Double-ended tube for Flip-Flop. Pack of 110 (100 + 10 spares).	


	Description	Flip-Flop Tube Seal (White Viton)	Part Number
			70-1304
	Details	White Viton tube seal for Flip-Flop (tube to funnel). Pack of 20.	


	Description	Flip-Flop Tube Seal (Perfluoroelastomer)	Part Number
			70-1306
	Details	Perfluoroelastomer tube seal for Flip-Flop (tube to funnel). Pack of 20.	


	Description	Funnel	Part Number
			70-2037
	Details	Funnel for Flip-Flop. Pack of 18.	


	Description	Funnel O-Ring	Part Number
			70-2039
	Details	O-ring for Flip-Flop funnel (funnel to vial cover). Pack of 20.	

	Description	Vial Seal (White Viton)	Part Number 70-1303
	Details	White Viton Vial Seal for Flip-Flop. Pack of 20.	

	Description	Vial Seal (Perfluoroelastomer)	Part Number 70-1305
	Details	Perfluoroelastomer Vial Seal for Flip-Flop. Pack of 20. NOTE: Notch in periphery indicates perfluoroelastomer.	


	Description	Ejector Peg	Part Number 70-1301
	Details	Vial Ejector Peg for Flip-Flop. Set of 18.	


	Description	Ejector Peg O-Ring	Part Number 70-1303
	Details	O-ring for Flip-Flop Vial Ejector Peg. Pack of 20.	

	Description	Vial Cover	Part Number 70-1300
	Details	Vial Cover for Flip-Flop. Set of 18.	


Dionex™ 250 ml ASE® Bottle Adaptor





	Description	Dionex 250 ml ASE Bottle Adaptor	Part Number
			70-2010
	Details	<p>Each bottle adaptor accepts one Dionex 250 ml ASE collection bottle (Part # 056284), allowing up to six bottles to be processed simultaneously. Set of six. Only compatible with Flask Rotors.</p>	


	Description	Dionex 250 ml ASE Bottle Adaptor Seal	Part Number
			70-2011/S
	Details	<p>Replacement seals for the Dionex 250 ml ASE Bottle Adaptor. Set of six.</p>	

Rocket Synergy - Accessories

	Description	Flask Rack	Part Number
	Details	Coated steel flask rack/carrier for safe storage or transport of six Standard Flasks, GC Vial System flasks or SampleGenie flasks.	
			70-1187


	Description	Flask Holder/Balance Nest	Part Number
	Details	Single flask stand for use as flask holder or balance nest. Compatible with Standard Flasks, GC Vial System flasks or SampleGenie flasks.	
			70-1797


	Description	Splash-Cap	Part Number
	Details	Viton splash-cap for use with Standard Flasks, GC Vial System flasks or SampleGenie flasks. Set of six.	
			70-1265

	Description	Seal Removal Tool	Part Number
	Details	Tool for safe removal of vial seal from SampleGenie or GC Vial System flasks without risk of damage to flasks or seals.	
			04-8321


Infinity Trolley





	Description	Infinity Trolley	Part Number
	Details	Infinity Trolley accommodates one Rocket Synergy evaporator, Julabo re-circulating chiller, and one five-litre waste container.	

	Description	Infinity Trolley Shelf Kit	Part Number
	Details	Shelf Kit for Infinity Trolley. Shelf accommodates one Edwards nXDS6i vacuum pump, or may be used as a sample preparation area.	


Rocket Synergy - General Spares


	Description	Outer Lid Seal	Part Number
	Details	Outer lid seal for Rocket Synergy, Rocket 4D Synergy and pre-Synergy Rocket evaporators.	
			04-5818/S


	Description	Sample Viewing Window	Part Number
	Details	Replacement sample viewing window and gasket assembly for Rocket evaporators with strobe sample viewing lights.	
			70-1182

	Description	Waste Solvent Container (5 Litre)	Part Number
	Details	Waste solvent container kit comprising five litre polypropylene waste container, two-port cap, drain hose (five metres) and fittings.	
			70-1247/S

	Description	Feed Lid Coupling O-Rings	Part Number
	Details	2 x Green perfluoroelastomer O-rings for feed lid joint coupling.	
			04-7943/S

	Description	Feed Hose	Part Number
	Details	70-1553/S Replacement 6 mm OD feed hose for Rocket 4D Synergy. 50 metre length.	

	Description	Hose Fittings	Part Number
	Details	70-1651/S Hose fittings kit comprising two Type 1A fittings and two Type 2N fittings for coolant, drain and exhaust hose connections.	

	Description	Hose-Tails	Part Number
	Details	70-2047/S Pair of hose-tails for connection to chilled water supply.	

Pre-Synergy Rockets - Rotor Compatibility

Early model Rocket Chemi evaporators (up to serial number RK0136) used a Flask Rotor with 'flask adaptors' attached. Some sample holders may not be compatible with Flask Rotors that still have flask adaptors attached.









Before using any sample holder for the first time, verify compatibility using the table on the following page.

If in doubt, please contact the Genevac Technical Support team at service@spscientific.com.



From left to right: Rocket Chemi rotor with flask adaptors, Rocket Chemi rotor with flask adaptors removed, and Rocket Synergy rotor.

The quick reference guide below shows the compatibility of sample holders with variants of the Rocket evaporator. **If in doubt, check the serial number of the evaporator to see if it is suitable for use with the sample holder.** Use of incompatible sample holders may result in damage to the evaporator.

Compatibility Quick Reference		Rocket Chemi (up to RK0136)	Rocket Chemi (RK0137 to RK0359)	Rocket Synergy (RK0400 onwards)
Sample holder type and variant				
Std Flask with Collar & Pad		<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Std Flask without Collar		<input checked="" type="checkbox"/>		
GC Vial System ²		<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
SampleGenie ⁴		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
SampleGenie ² SampleGenie ³ (discontinued)		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
SampleGenie ¹ (discontinued)		<input checked="" type="checkbox"/>		
Puck/Flip-Flop		<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Dionex 250 ml ASE Bottle Adaptor				<input checked="" type="checkbox"/>

<input checked="" type="checkbox"/>	Compatible
<input checked="" type="checkbox"/>	Low Rotor Speed upgrade is recommended.
<input type="checkbox"/>	Sample holder may only be used if flask-adaptors have been removed from the rotor.

Read these instructions carefully and keep them near the system for ease of reference. Your attention is drawn, in particular, to the **Limitations of Use** and **Rotor Compatibility** sections.

These instructions are correct at the time of going to press and may be subject to change without notice. Some of the features described within this manual may not apply to equipment manufactured before this manual's publication date; this includes systems that have been upgraded.

No part of these instructions may be reproduced in any form or be processed, duplicated or distributed by electronic or other means without the express written permission of Genevac Limited.

Should you need to contact Genevac for assistance, please use one of the contact methods shown. Please have the equipment serial number and any other pertinent information to hand.

Sample holders should not be discarded in your regular waste stream. Contact your representative or Genevac for proper disposal instructions.

Within the EU, it is Genevac's responsibility under the WEEE Directive to provide for the recycling of Genevac products.

Genevac Ltd.
The Sovereign Centre
Farthing Road
Ipswich
IP1 5AP
United Kingdom

Genevac Inc.
SP Scientific
3538 Main Street
Stone Ridge
NY12484
United States of America

Tel:

Sales: +44 (0) 1473 240000
Service: +44 (0) 1473 243000
Fax: +44 (0) 1473 461176

Email:

Sales: salesinfo@genevac.co.uk
Service: service@spscientific.com

Application Support: applications@genevac.co.uk

<http://www.spscientific.com>

Tel:

Sales: (1) 845 267 2211
Service: (1) 877 548 4666
Fax: (1) 845 267 2212

Email:

Sales: sales@genevacusa.com
Service: service@spscientific.com

Application Support: applications@genevac.co.uk

<http://www.spscientific.com>

